

Qualy Bone Medical



75%

HYDROXYAPATITE

Z5%
TRICALCIUM

PHOSPHATE (B-TCP)

- POROUS SYNTHETIC
 CERAMIC DESIGNED
 FOR THE FILLING OF
 BONE VOIDS OR
 DEFECTS.
- REGENERATION
 AND BONE GROWTH,
 STIMULATING THE
 PROLIFERATION AND
 DIFFERENTIATION OF
- DUE TO ITS MESH
 CONFIGURATION IT
 HAS A HIGH POROSITY
 FACILITATING THE
 PROLIFERATION OF THE
 BONE CELLS, INTO THE
 VOIDS
- THIS FORMAT GIVES IT EXCELLENT MECHANICAL RESISTANCE
- CELL ADHESION
 OBSERVED
 AFTER 5 DAYS
 OF BONE
 GRAFTING

DICATIONS

- FILLING DEFECTS OR BONE CAVITIES
- VOLUME INCREASE IN BONE DEFECTS THAT ARE NOT INTRINSIC TO THE STABILITY OF THE BONE STRUCTURE

Fractures with bone defect Vertebral arthrodesis Tibial osteotomy Tibial and femoral fracture Total knee and hip revision Spine Surgery



E QUALYBONE BO

SAFETY

100% SYNTHETIC AND RESORBABLE

HIGH COHESIVITY AND VOLUME CONSERVATION

BEING VERY HYDROPHILIC GIVES IT AN EXCELLENT COHESIVENESS OF THE PARTICLES, IN TURN ALLOWS THE CONSERVATION OF THE INITIAL CAVITY **VOLUME**

RESORBABLE

FROM 6 TO 24 MONTHS, DURING THIS PERIOD IT IS REPLACED BY NEW BONE TISSUE

VASCULARIZATION


INDUCES A REMARKABLE VASCULARIZATION

RADIOPAQUE

EXCELLENT PROPERTY FOR RADIOLOGICAL CONTROL DURING THE PROCESS OF BONE REGENERATION UNTIL ITS TOTAL RESORPTION

DIRECT APLICATION

NO MEMBRANE REQUIRED. DUE TO ITS PHYSICAL AND MECHANICAL PROPERTIES



GEOMETRY	SIZE	QUANTITIES	REF.
	3 - 4 MM	2,5 G	QBB030402G
	3 - 4 MM	5 G	QBB030405G
	3 - 4 MM	10 G	QBB030410G
	3 - 4 MM	15 G	QBB030415G
GRANULES	3 - 4 MM	20 G	0BB030420G
	3 - 4 MM	25 G	0BB030425G
	3 - 4 MM	30 G	0BB030430G
	4 - 7 MM	2,5 G	0BB040702G
	4 - 7 MM	5 G	0BB040705G
	4 - 7 MM	10 G	0BB040710G
	4 - 7 MM	15 G	0BB040715G
	4 – 7MM	20 G	QBB040720G
	4 – 7MM	25 G	QBB040725G
	4 – 7MM	30 G	QBB040730G
BLOCKS	4X4X20 MM	1 UNIT	QBB040420B
	8X8X20 MM	1 UNIT	QBB080820B
	15X15X20 MM	1 UNIT	QBB151520B
	15X20X30 MM	1 UNIT	QBB152030B
CYLINDER	8X8X20 MM	1 UNIT	QBB080820C
WEDGE	6X25X30 MM 8X25X30 MM 10X25X30 MM 12X25X30 MM 14X25X30 MM	1 UNIT 1 UNIT 1 UNIT 1 UNIT 1 UNIT 1 UNIT	QBB062530W QBB082530W QBB102530W QBB122530W QBB142530W
STICK	4X4X20 MM	4 UNITS	QBB040420Q
	4X4X20 MM	5 UNITS	QBB040420F



Qualy Live

Rua Dona Ines de Castro, 17C - Bloco B sala 15

2720-193 • Amadora • Portugal

⊕ +351 215861750

+351 935652996

 □ qualylive@qualylive.com www.qualylive.com

